

Determination of Vitamin D in Vitamin D Powder

1. Principle

The microencapsulated vitamin D powder, after being diluted and broken, is ultrasonically extracted with absolute ethanol and then determined by high performance liquid chromatography.

The vitamin D reference standard is used as a control and the external standard method is used for quantification.

2. Reagent

2.1 0.01 mol/L hydrochloric acid solution: Draw 0.86 ml of hydrochloric acid solution(AR), and add water to 1000 ml.

2.2 Anhydrous ethanol: Analytical reagent (deacetalized)

2.3 Methanol: HPLC grade

2.4 Ultra-pure water

2.5 Standard stock solution:

Accurately weigh 0.02g (accurate to 0.0001g) of the reference standard by purity, equivalent to 40,000,000IU/g vitamin D.

Dissolve in absolute ethanol and dilute to a 50-ml brown volumetric flask and shake well.

2.6 Standard solution:

Accurately draw standard stock solution 5.00ml in a 250-ml brown volumetric flask, dilute with absolute ethanol to volume.

Shake well, pass the solution through an organic syringe filter of 0.22- μ m pore size, and take the filtrate for injection test. (Available at the time of use).

3. Sample processing

Accurately weigh 0.4g (accurate to 0.0001g) sample in a 250 ml brown volumetric flask and add 10ml of 0.01mol/L hydrochloric acid solution, treated in an ultrasonic bath at 60 ° C for 15 min to completely disintegrate, and then add anhydrous ethanol to the volume.

Shake well, place it in an Ultrasonic oscillator for 10 min, remove, cool to room

temperature, dilute with absolute ethanol to volume.

Shake well, pass successively the solution through organic syringe filter of 0.45µm and 0.22-µm pore size, and take the filtrate for injection test.

4. Chromatographic conditions

Column: C18 column, 250mm × 4.6mm, particle size: 5µm

Mobile phase: methanol

Flow rate: 1 ml/min

Detection wavelength: 254 nm

Injection volume: 20µL

5. Determination of previtamin D correction factor

Accurately weigh the standard stock solution 5ml to 100ml saponification flask, add 3 BHT and appropriate amount of absolute ethanol, and reflux in a water bath at 90 ° C for 45 minutes, then cool down. Transfer to a 50-ml brown volumetric flask with an appropriate amount of absolute ethanol, add absolute ethanol to volume and shake well.

Pass through an Organic filter membrane of 0.22- µ m pore size, and use high-performance liquid chromatograph as the above-mentioned chromatographic conditions for injection detection to obtain the the peak area of vitamin D **A'** and the peak area of pre-vitamin D **A'pre**.

The standard solution is also passed through an Organic filter membrane of 0.22- µ m pore size, and use high-performance liquid chromatograph as the above-mentioned chromatographic conditions for injection detection to obtain the peak area of vitamin D **A** and the peak area of pre-vitamin D **Apre** .

Calculation formula:

Calculate the Pre-vitamin D correction factor (F):

$$F = (A - A') / (A'_{pre} - A_{pre})$$

6. Calculation formula:

$$X \text{ (IU/g)} = \frac{(A_0 \times F + A_1) \times M_2}{A_2 \times M_1 \times 10} \times 40000000$$



HEALTH SOURCES
NUTRITION CO.,LTD.

X—the content(IU/g) of Vitamin D in the sample

A 0—the peak area of Pre-vitamin D in the sample

A 1—the peak area of Vitamin D in the sample

A 2—the peak area of Vitamin D in the standard

M 1—Sample weight

M 2—Standard weight

7. Precision

The absolute difference between two independent determinations obtained under repeatability conditions shall not exceed 5% of the arithmetic mean.