

# Heat-labile UDG

Catalog # P051-01



Version 5.1

Vazyme biotech co., ltd.

## Introduction

UDG (Uracil-DNA Glycosylase) catalyses and hydrolyses uracil base and N- glycosidic bond of the sugar phosphate skeleton which belong to single or double-stranded DNA containing dU. Thus, free uracils are released and the base-free sites are dissociated easily by hydrolysis. Heat-labile UDG stemming from psychrophilic marine bacterium, is sensitive to high temperature and therefore inactivates enzymes irreversibly when temperature over 50°C, which is suitable for PCR/QPCR, RT- PCR /RT-QPCR system.

## Package Information

Components	P051-01 100U
Heat-labile UDG (1 U/μl)	100 μl

## Storage Buffer

20 mM Tris-HCl, pH 8.0@ 25°C  
0.1 mM EDTA  
100 mM KCl  
1 mM DTT  
50% Glycerol (v/v)  
0.5% NP-40 (v/v)  
0.5% Tween-20 (v/v)

## Storage

Store at -20 °C

## Origin

Recombinant *E. coli* strain cloned with UDG gene of psychrophilic marine bacterium

## Unit Definition

One One unit (U) is defined as the amount of enzyme that releases 1 nmol of uracils from the DNA strand (containing dU) within 1 hour at 37°C in the reaction system containing 70 mM of Tris - HCl, pH 7.5, 10 mM of NaCl, 1 mM of EDTA, 100 μg/ml of BSA reaction liquid.

## Specific Activity of Enzyme

≥200,000 U/mg

## Quality Control

Exonuclease residue detection: DNA electrophoresis bands do not change when 10 U of this enzyme and 0.6 μg of λ-Hind III are incubated at 37°C for 16 hours.

Endonuclease residue detection: DNA electrophoresis bands do not change when 10 U of this enzyme and 0.6 μg of Supercoiled pBR322 DNA are incubated at 37°C for 4 hours.

RNase residue detection: RNA electrophoresis bands do not change when 10 U of this enzyme and 1 μg of total RNA of HeLa cell are incubated at 37°C for 1 hours.

*E. coli* DNA residue detection: *E. coli* genome residue of 200 U of this product should be less than 10 copies in TaqMan qPCR detection specified with *E. coli* 16s rDNA .



Vazyme Biotech Co., Ltd  
www.vazyme.com

Order: global@vazyme.com

Support: global@vazyme.com

**For research use only, not for use in diagnostic procedures.**

## Protocol

### 1. Recommended reaction mixture for PCR:

ddH <sub>2</sub> O	to 50 µl
10 × Taq Buffer( with 20 mM MgCl <sub>2</sub> )	5 µl
25 mM MgCl <sub>2</sub> <sup>a</sup>	Optional
dUTP <sup>b</sup>	0.6 mM
dATP/dCTP/dGTP	0.2 mM each
Template DNA	Optional
Primer 1 (10 µM)	2 µl
Primer 2 (10 µM)	2 µl
Taq DNA Polymerase (5 U/µl)	0.5µl
Heat-labile UDG (1 U/µl) <sup>c</sup>	1 µl

a. The final concentration of Mg<sup>2+</sup> can be adjusted between 2.0 and 3.0 mM according to experiment needs.

b. The final dUTP concentration can be adjusted to 0.2 - 0.6 mM according to experiment demands.

c. According to the requirements of the experiment, the general amount of the Heat-labile UDG in the 50 µl reaction system is 0.1 - 1 U.

### 2. PCR conditions:

25 °C	10 min	U-containing template degradation
95 °C	2 min	UDG inactivation, template degeneration
PCR Reaction		

