Wayen Exosome Isolation Kit

User Manual

Cat# EIQ3-03001 (Urine) Version 2017-01

Description (For research only)

EIQ3-03001 kit is used to isolate exosomes between 30 and 200 nm diameter from urine. By adding appropriate amount of reagents (A and B) to urine sample, and incubating the mixture over 12 h, the exosomes can be collected with centrifugation.

Advantages

- ✓ Quantity: Higher yield (versus other kits or methods)
- ✓ Easy: No ultracentrifugation

Contents

EIQ3-03001 kit contains Reagent A and B.

Storage

Reagent A is shipped at room temperature and Reagent B is shipped at 4 °C. The kits should be stored properly after receipt. Properly stored kits are stable for 1 year from the date received.

Components	Storage	Volume
Reagent A	Room Temperature	37.5 mL
Reagent B	-20 °C	3.35 mL

Experiment Protocol of EIQ3-03001 (Urine)

1. Prepare Urine Sample

1.1 Take the urine sample from storage and keep it on ice. If starting with frozen sample, thaw the sample completely in a 37 °C water bath and then place it on ice.

1.2 Centrifuge the urine sample at $3,000 \times g$ for 15 minutes at 4 °C.

1.3 Transfer the supernatant to a new tube and place it on ice until ready to perform the isolation.

2. Isolate Exosomes (Balance the Reagent A to room temperature before use and the starting volume of urine is recommended to be 20 mL. The protocol below is shown with 20 mL urine.)

2.1 Take out 20 mL pre-treated urine, add 7.5 mL Reagent A, mix the mixture well by inverting it up and down until obtain a homogenous mixture.

Note: For every 10 mL of urine, 3.75 mL Reagent A should be added. When the start urine volume is changed, the volume of the reagent A should also be changed.

2.2 Add 670 μ L Reagent B to the mixture above and invert the tube up and down to obtain a homogenous mixture.

Note: For every 10 mL of urine, 335 μ L Reagent B should be added. When the start urine volume is changed, the volume of the reagent B should also be changed.

2.3 Incubate the mixture at 4 °C overnight (12-16 h).

2.4 After incubation, centrifuge at $3,000 \times g$ for 60 minutes at 4 °C.

2.5 Take out 1 ml supernatant into a 1.5 mL tube firstly and then remove the residual supernatant.

2.6 Resuspend the pellet completely with the 1 ml supernatant above. Transfer the mixture to a new 1.5 mL tube.

2.7 Centrifuge the re-suspension at $10,000 \times g$ for 10 minutes at 4 °C. Remove the supernatant without disturbing the precipitated pellet.

2.8 Resuspend the pellet with 50-200 μ L $1 \times$ PBS, and mix it well by vortexing or pipetting up and down until obtain a homogenous mixture.

2.9 Centrifuge the sample at $10,000 \times g$ again for 5 minutes at 4 °C. Transfer the supernatant to a new tube.

2.10 The supernatant contains exosomes. The exosomes can be used for downstream analysis immediately or aliquoted and stored at -80 °C till next experiment.

Notice

This kit is for research use only, not for clinical diagnostic purpose.

$1 \times$ sterile PBS is not supplied and should be prepared by user.

We recommend that exosomes used for electron microscopy, NTA analysis and proteomics studies should be filtered by 0.22 μ m filtration.

For more detail information please visit our official website: www.wayenbio.com

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