




























1.1 Symbols Explained

| Symbol | Explanation | Symbol | Explanation |
|---|--|--|--|
|  | Type B applied part |  | OFF (power, Scanner) |
|  | Alternating current |  | ON (power, Scanner) |
|  | The equipment should be disposed of according to local regulation for disposal of electronic equipments. |  | Class 1 laser product |
|  | Protective earth (ground) |  | Refer to instruction manual/booklet |
|  | Catalogue number |  | Serial number |
|  | Manufacturer |  | Authorized representative in the European Community |
|  | Date of manufacture |  | Medical Device |
|  | Use-by date 使用期限 |  | Unique device Identifier |
|  | Crushing of hands | | European union approval |
|  | This Way Up |  | “OFF” for part of the equipment |
|  | Fragile, handle with care |  | “ON” for part of the equipment |
|  | DO NOT STACK |  | DO NOT ROLL |
|  | Keep Dry |  | Temperature Limit |
|  | Humidity Limitation |  | Atmospheric Pressure Limitation |
| Rx ONLY | Caution: Federal law prohibits dispensing without prescription. | voltage selector | The selector switch of input voltage, which can be toggled to adapt to the network supply voltage. |