

PDL1/aAPC/Raji

CBP74145

Contents

I. Background.....	1
II. Introduction.....	1
III. Representative Data.....	2



PDL1/aAPC/Raji CBP74145

I. Background

N/A

II. Introduction

Expressed gene: PD-L1、 aAPC

Stability: 32 passages (in-house test, that not means the cell line will be instable beyond the passages we tested.)

Freeze Medium: 90% FBS+10% DMSO

Culture Medium: RPMI-1640+10%FBS+1ug/ml puromycin+200ug/ml hygromycin

Mycoplasma Testing: Negative

Storage: Liquid nitrogen

Application(s): Functional(Report Gene) Assay



III. Representative Data

Dose Response of Anti-PD-L1 Blocking Antibody in PD-1-IL2 Pathway Effector Reporter Cells (C37) With Raji PD-L1 aAPC Cells




Figure 1. Dose Response of Anti-PD-L1 Blocking Antibody in PD-1-IL2 Pathway Effector Reporter Cells (C37) With Raji PD-L1 aAPC Cells.

Dose Response of Anti-PD-1 Blocking Antibody in PD-1-IL2 Pathway Effector Reporter Cells (C37) With Raji PD-L1 aAPC Cells




Figure 2. Dose Response of Anti-PD-1 Blocking Antibody in PD-1-IL2 Pathway Effector Reporter Cells (C37) With Raji PD-L1 aAPC Cells.

